


## Adjustable Socket Order Form

<b>Bill To:</b> Address: _____ _____ _____  <b>Ship To:</b> Address: _____ _____ _____  <i>Same as Bill To</i>	<b>Patient Name:</b> _____  Height: _____ <input type="checkbox"/> Left <input type="checkbox"/> Male Weight: _____ <input type="checkbox"/> Right <input type="checkbox"/> Female Age: _____ <input type="checkbox"/> Bilateral  K Level:    1        2        3        4        Leg    Hip  <b>Order Date:</b> _____  <b>In-Office Request Date:</b> _____ <input type="checkbox"/> am <input type="checkbox"/> pm  <b>PO / OPS #:</b> _____
<b>Practitioner:</b> _____ Phone #: _____	


## Socket Configuration


Please indicate the type of socket to be fabricated.





**Transfemoral Construction Details**  
**Standard:** Proximal strap for volume adjustment  
 Option: Add a second zone that includes adjustable panels (choice C). Discuss with technician

**Select reel location:**  
**Standard:** A-Proximal mount  
 Options: B-Distal mount or C-add second zone; discuss with technician

  
A Transfemoral - Standard

  
B Transfemoral - Distal reel

  
C Transfemoral multi-zone



**Transtibial Construction Details**  
**Standard:** 3 moveable panels: 2 anterior & 1 posterior  
 Options: Select the number of panels you prefer

**Select panel location (provide size and shape) and reel location:**

**Panel Choice(s)**


Anterior Lateral   
  Anterior Medial   
  Posterior

**Options:**

- One posterior panel for muscle expansion of the calf.
- Two anterior panels for creating rotation control.
- Three panels, posterior panel eases donning/doffing; anterior panels manage volume.

**Reel Location**

Socket / Frame   
  Panel





**Symes Construction details:**  
**Standard:** 1 moveable panel  
 Options: locate laterally, medially, or posterior

**Select panel location (provide size and shape) and reel location**

**Panel Choice**

Lateral Panel   
  Medial Panel   
  Posterior Panel

**Reel Location**

Socket / Frame   
  Panel