

# INTERSECT

## Custom Leather AFOs

High-Tech & Comfort Converge



**NATIONAL LABS**  
ORTHOTICS • PROSTHETICS • INNOVATION

# INTERSECT. Custom Leather AFOs.

High Tech & Comfort Converge.

**Advanced Laser-scan modeling** provides

- Custom-tailored accuracy
- Digital historical patient data retention
- Fast track ordering & -2 week standard turnaround time

**Hand-crafted & wide selection** of options to optimize each patient's unique clinical goals and provide supreme comfort

## SIX BASIC DESIGNS

### STANDARD

- 8" height

### EXTENDED

- Standard with proximal extension

### SYNTHETIC

- Manmade material with polycore

### RIGID

- Reinforced with solastic, no liner

### SOFT

- Plastazote & aliplast with poly core

### ARTICULATED

- Any of the above designs with articulating ankle (Tamarack standard)

## OPTIONS

- Full-length foot plate
- Velcro, speed laces or boot hooks combinations.
- Plastic reinforcements can be integrated for SMO trim or solid ankle designs.
- Partial foot configuration is available.
- Colors: sand, black, brown or white
- Leather (standard) or
- Synthetic
- BOA Closure